

# QuadLink™ Presutured Construct

## Product Highlights

The QuadLink construct is a presized and preassembled allograft quadriceps tendon presutured with #2 FiberLoop with FiberTag™ suture, designed to be used with ACL TightRope® implants.

### Features and Benefits

- Minimal graft preparation time
- Presized diameter and lengths
- Sterile via Allowash XG®\* technology
- Preassembled with #2 FiberLoop with FiberTag suture
- Preloaded with passing sutures to facilitate loading with ACL TightRope II implants

The all-inside ACL technique is ideal for quadriceps grafts, as only a graft length of 60 mm to 75 mm is needed. ACL TightRope II implants are used to fixate the graft in minimally invasive sockets created with a FlipCutter® III drill. A standard all-inside technique is used to implant the graft.




## Ordering Information

Diameters available: 9.0 mm to 11.0 mm (in .5 mm increments). Lengths available: 60 mm to 75 mm (in 1 mm increments). Tendons ordered directly from LifeNet Health® (1-888-847-7831).

Product Description	Item Number
<b>Implants</b>	
QuadLink™ Presutured Construct	<b>FQL</b>
TightRope® II BTB, Recon IB™	<b>AR-1588BTB-IB</b>
TightRope II RT, w/ double-loading passing suture	<b>AR-1588BTB-2J</b>
TightRope II ABS, implant open	<b>AR-1588TN-21</b>
Concave ABS Button, 11 mm w/ 4 mm collar	<b>AR-1588TB-3</b>
Concave ABS Button, 14 mm w/ 7 mm collar	<b>AR-1588TB-4</b>
Concave ABS Button, 20 mm w/ 9 mm collar	<b>AR-1588TB-5</b>
<b>ACL Backup Kits</b>	
Implant System, secondary fixation w/ BioComposite SwiveLock® anchor, 4.75 × 19.1 mm	<b>AR-1593-BC</b>
Implant System, secondary fixation w/ PEEK SwiveLock Anchor, 4.75 × 19.1 mm	<b>AR-1593-P</b>
<b>Instruments for FlipCutter® Technique</b>	
RetroConstruction™ Drill Guide Set	<b>AR-1510S</b>
Footprint Femoral ACL Guide, right	<b>AR-1510FR</b>
Footprint Femoral ACL Guide, left	<b>AR-1510FL</b>
FlipCutter III Drill	<b>AR-1204FF</b>
<b>For Medial Portal Technique</b>	
Transportal ACL Guides (TPGs), 4 mm-8 mm	<b>AR-1800-04 – 08</b>
Low-Profile Reamers, 5 mm-13 mm (in .5 mm increments)	<b>AR-1405LP – 1413LP</b>
TightRope Drill Pin, open	<b>AR-1595T</b>
TightRope Drill Pin, closed	<b>AR-1595TC</b>

Product Description	Item Number
<b>Accessories</b>	
Suture Retriever	<b>AR-12540</b>
Suture Cutter for ACL TightRope Suture	<b>AR-4520</b>
Graft Sizing Block	<b>AR-1886</b>
Graft Tube Set	<b>AR-1886S</b>
<b>GraftPro® Graft Preparation System (AR-2950DS) Includes:</b>	
GraftPro Board	<b>AR-2950D</b>
GraftPro Posts	<b>AR-2950AP</b>
GraftPro Case	<b>AR-2950DC</b>
GraftPro GraftLink® Tensioner	<b>AR-2950GT</b>
GraftPro GraftLink Holder	<b>AR-2950GH</b>
GraftPro Button Holder	<b>AR-2950BH</b>
GraftPro Soft-Tissue Clamp	<b>AR-2950SC</b>
<b>Suture</b>	
FiberStick™ Suture, #2 FiberWire® Suture, 50 in (blue), one end stiffened	<b>AR-7209</b>
TigerStick® Suture, #2 TigerWire® Suture, 50 in (white/black), one end stiffened	<b>AR-7209T</b>
0 FiberLink Suture, FiberWire Suture w/ 1.5 in closed loop at one end (blue)	<b>AR-7258</b>

## Loading of the femoral graft end with TightRope BTB implant



This description of technique is provided as an educational tool and clinical aid to assist properly licensed medical professionals in the usage of specific Arthrex products. As part of this professional usage, the medical professional must use their professional judgment in making any final determinations in product usage and technique. In doing so, the medical professional should rely on their own training and experience and should conduct a thorough review of pertinent medical literature and the product's directions for use. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes.

Products may not be available in all markets because product availability is subject to the regulatory approvals and medical practices in individual markets. Please contact Arthrex if you have questions about the availability of products in your area.

View U.S. patent information at [www.arthrex.com/corporate/virtual-patent-marking](http://www.arthrex.com/corporate/virtual-patent-marking)  
 \*Allowash XG is a registered trademark of LifeNet Health.

**arthrex.com**

© 2021 Arthrex, Inc. All rights reserved. FL1-000037-en-US\_C

