



GraftNet™ XL Bone Collection Device

Product Highlight


The GraftNet XL device is designed to maximize autologous bone graft collection and uses a large-volume mesh collection filter to easily collect autologous bone to integrate the patient's own bone graft into procedures.

Features and Benefits

- 1. Collect bone in line with desired shaver and suction:** Universal adaptors at inflow and outflow ensure the GraftNet XL device can be quickly incorporated into any procedure.
- 2. Quick access to high-volume bone graft:** Ergonomically designed lid with a lock-and-unlock feature to protect and provide access to the collected bone graft.
- 3. Easily integrate into sterile field workflow:** The drape clip allows the GraftNet XL device to be stably fixed to the sterile field. The assembled tubing clamp is available to control suction flow as needed throughout the procedure.



As Simple as Resect and Collect™: GraftNet XL With the BioACL™ Technique



1 Attach GraftNet XL tubing to the desired suction.



2 Collected bone is filtered through the internal filter of the GraftNet XL device.



3 Once bone collection is complete, disconnect the device from suction and move the cap from the locked to the unlocked position.



4 Remove and invert the filter to free collected bone into a sterile basin. A Freer elevator can be used to aid in removal of tissue from the filter.

Ordering Information

Product description	Item number
GraftNet XL bone collection device	ABS-10052

arthrex.com

© 2024-11 Arthrex, Inc. All rights reserved. FL1-000806-en-US_B

