

Mini CFS

Mini Comprehensive Fixation System

1.4 mm/1.6 mm/2.0 mm/2.4 mm




1.4 mm




1.6 mm



2.0 mm



2.4 mm



Mini Plates...Big Solutions

- Comprehensive solution for mini fragment fixation and fusion
- 56 plating options
- Color-coded modules with corresponding instrumentation
- Radiolucent and optimized instrumentation

Arthrex®



Mini CFS

Comprehensive Plating Options




Z-Plate
1.4 mm/1.6 mm



Hook Plate
1.4 mm/1.6 mm



Cage Plate
1.4 mm/1.6 mm/2.0 mm/2.4 mm



Bridge Plate
1.4 mm/1.6 mm/2.0 mm/2.4 mm



Straight Plate
1.4 mm/1.6 mm/2.0 mm/2.4 mm



Y-Plate
1.4 mm/1.6 mm/2.0 mm/2.4 mm



T-Plate
1.4 mm/1.6 mm/2.0 mm/2.4 mm




Triangle Plate
1.4 mm/1.6 mm/2.0 mm/2.4 mm



Rotation Correction Plate
1.6 mm/2.0 mm



Straight Plate, Reinforced
2.0 mm/2.4 mm



V-Plate
2.0 mm/2.4 mm



Wide T-Plate
2.0 mm/2.4 mm



Corresponding Screw Fixation




Low Profile Screw, Cortical, Titanium
1.0 mm*: 6 mm-14 mm
1.4 mm: 6 mm-16 mm
1.6 mm: 6 mm-24 mm
2.0 mm/2.4 mm: 6 mm-40 mm



Variable-Angle Locking Screw
1.4 mm: 6 mm-16 mm
1.6 mm: 6 mm-24 mm
2.0 mm/2.4 mm: 6 mm-40 mm



Micro Compression FT™ Screw
2.5 mm – Optional caddy for intramedullary fixation
36 mm-50 mm



Mini Compression FT™ Screw
3.5 mm – Optional caddy for intramedullary fixation

*FDA clearance pending

www.arthrex.com

© 2020 Arthrex, Inc. All rights reserved. FL1-000163-en-US_A